

Tethered Lean & Release System MATLAB GUI Instructions and Design

Notes - LSL based

Tuesday, October 29, 2024 10:35 AM

App Design

Purpose: To display and plot the tether tension/estimated lean torque/estimated lean angle in real time. To provide researcher a clear indication of whether the test subject has lean enough/ready to be released.

Functionalities:

Calculate lean angle and lean torque based on the force sensor readings and input parameters.

Plot results in real time.

- Initiate: Initialize the library, obtain available channels based on channel properties
- Start: connect to the channel inlet and start receiving data/time stamp vectors
- Stop: stop streaming and terminates the initialization
- Info:**
 - a. Stream Info: Prop type: the type of the stream property, such as "name", "type", "source_id", etc
 - b. Prop name: the name of the stream property, such as "Name" -> "MyKeyboard"
 - c. Pull type: stream or chunk
 - d. Estimated type: lean angle or lean torque
 - e. Target Value: toggle between target lean angle and target lean torque, depending on the estimation type
- Weight: a. Weight b. height
- Auxiliary info: a. height-CoM coefficient: the height of the CoM of patient over the height of the patient ratio, choose between established average or user-defined
- b. tether height
- Calibrate: calibrate the force sensor with known weights
 - a. Calib Weight: input the known weight
- Log: display detailed info and error messages


Save File: save recorded data.

Earlier Versions Interface Layout



Use Instructions:

- Preparation:
 - Place the subject on a harness.
 - Connect the harness to BrainVision Recorder, etc.
 - Put the harness on the subject.
 - Describe the experiment instructions to the subject.
 - Record the subject's height in cm.
 - Record the subject's weight in kg.
 - Note the subject's sex.
- Setup the tethered-release system:
 - Stretch the harness to about 250cm over the height of the metal clip on the subject's harness when the subject is standing straight with the harness on. Record the pulse height in cm.
 - Clip the tether to the bice release and the clip on the harness.
 - Ask the tethered subject to walk away from the tethered-release system, following the red tape on the ground, until the tether is stretched out.
 - According to the experiment requirements, ask the subject to step back 1 or 2 steps to learn.
- Setup the LiveAmp:
 - Switch on the BrainVision LSL Connector on the Desktop of the Moll Laptop.
 - Linkamp Connector
- Press Scan. The Target LiveAmp should appear in the Available Devices dropdown. We can confirm the LiveAmp by comparing the number of EEG channels will automatically change to 64 if we are using 2 LiveAmps.
- Specify the LSL Channel for our STE box, we type in 8.
- Check the Tetherable ADC connection.
- Keep the Chest and Sampling Rate default. → 500Hz R = 100Hz raw data.
- Check the Channel under LSL Trigger Output Style. This allows a standalone EEG channel-style marker for de-jittering in later analysis. (However, this functionality and its compatibility to the MATLAB App hasn't been fully checked as of 04/22)
- If successfully linked, the Link button will change to Unlink.
- Setup the Lab Recorder:
 - Open the Lab Recorder. The Lab Recorder is used to record data, but we can also check important information of the stream.



- If the LiveAmp and TriggerBox and STE are already streamed. We should be able to see the streams on the left side Record from Streams box. The green text inside the parentheses corresponds to the name of the stream. In the case of LiveAmp, TriggerBox, and STE connected, should be able to see 3 streams. The name of the stream will be the device followed by DeviceTrigger or STETrigger.
- Click Select All.
- If changes made to the streams, click Update. The updated streams will show up in the Record from Streams box shortly.
- Click Save path and save file name.
- Click Start to start recording. But you may want to start recording after using the MATLAB App to visualize all important sensors and triggers and make sure that they function properly.

- Setup the MATLAB App (Tethered_Lean_Release_Feedback):
 - On the MATLAB, open MATLAB R2024a. Select the APPS tab, the Tethered_Lean_Release_Feedback app is the 1st app to the right of the Package App icon.

b. Click the Tethered_Lean_Release_Feedback app to switch on the app. You will see the following panel, the same as the layout shown above.

- The marker stream will automatically find the type marker stream. Specify the prop type and value similarly.
- The log will display the founded stream. Cross validate the stream information to make sure that we connected to the correct stream.
- Configure Stream Info: the app can subscribe to 1 sensor and 1 marker stream. Depending on your need, specify the Sensor Stream Info and Marker Stream Info to find specific streams.
- We can select Prop type to find streams according to their specific properties, such as their name, type, or source_id.
- For instance, we can find the LiveAmp EEG Stream by: (check the name of the streams via the Lab Recorder)

- Sensor Prop Type: EEG
- Sensor Prop Value: (054207-156<-->)

→ QDO : How to get the marker stream.

- The marker stream will automatically find the type marker stream. Specify the prop type and value similarly.
- The log will display the founded stream. Cross validate the stream information to make sure that we connected to the correct stream.

d. Configure Stream Info:

- Input the patient weight and height.

i. Patient Info:

- Weight (kg) [60]
- Height (cm) [170]

ii. Configure Aux Info:

- Input the Tether Height (which is the metal plate height from the ground) in the previous "b. Setup the tethered-release system." stage.

iii. Configure the CoM type according to the subject's sex. You can specify your own CoM ratio by selecting custom in the CoM Type and type in the CoM.

iv. Com Coefficient: Male Avg → MALE .65 FEMALE .55 → are legal!

v. Configure Stream Controls:

- Switch off Marker Switch if you do not have/want to stream a marker. If no marker is detected, the marker switch will automatically turn off.
- Marker Switch: Off [] On
- Initiate button will check for streams that satisfy the prop type and values specified in the "c. Stream Info" section.

vi. See Log for details of the connected stream or if the stream is connected.

Log:


```
>Loading the library
Sensor Stream Info
Marker Stream Info
Ready to Release
Sensor Prop Type: MyAudioStream
Marker Prop Type: Keyboard
Sensor Prop Value: MyAudioStream
Marker Prop Value: Keyboard
Patient Info
Auxiliary Info
Calibrate
Weight (kg): 60
Calib Weight (kg): 0
Height (cm): 170
Calib Type: Male Avg
Calib Coefficient: 0.5
Tether Height (cm): 150
Stream Controls
Initiate
Start
Terminate
Full Type: chunk
Marker Switch: Off [ ]
On
Log
File Name: Recording_sample_jyyy_MM
File Type: mat
Save Path: Working Folder
Save
```

vii. Start Button will start the plotting/recording. The sensor values are shown on the right. You can select a specific channel to plot by inputting different channel numbers in the Channel Selection.


Log:


```
MATLAB App
Sensor Stream Info
Marker Stream Info
Ready to Release
Raw Readings Type: mV
Raw Readings: 40
Raw Voltage Readings
Estimated Lean Angle
```

viii. Start Button will start the plotting/recording. The sensor values are shown on the right. You can select a specific channel to plot by inputting different channel numbers in the Channel Selection.

Log:


```
MATLAB App
Sensor Stream Info
Marker Stream Info
Ready to Release
Raw Readings Type: mV
Raw Readings: 1428
Raw Voltage Readings
```



Version History

Updated 12/10/2023

- Initiate -> can initiate library engine.
- Stop -> can stop streaming and delete stream info anytime.
- Start -> can successfully create inlet and receive data.
- Pull type toggle between chunk and sample, should use chunk by default to avoid missing data points
- Log is correctly displayed
- Raw data is correctly displayed
- Estimated data is correctly plotted
- Estimated type is correctly toggled
- Scrolling display

- Patients info -> weight and height is correctly recorded and used
- Aux Info CoM coefficient is correctly selected and used
- Calibration is correct -> needs to build a rack for the known weight to direct tension to gravity direction, see "Force Sensor Zero and Span Calibration Setup"

- Channel selection -> with a multi-channel input, find the correct channel with data, need to be tested with the EEG cart ActiChamp
- Indicator status determination -> based on avg data values over a time period instead? To ensure stability
- Display target lean angle and target lean torque threshold?
- Calibration
- Save data
- Test for performance, latency and consistency, etc.

Update 12/12/2023

Timers object causes crashes, consider implementing with the newer time table object instead.

To Do:

- Implement time table
- Object declaration/initialization
- Object append new data
- Append new data
- Test for performance, latency and consistency, etc.

Update 04/21/2024

1. Check the voltage to force conversion



2. Check the force to angle conversion

3. Swapped Controls and Info Panel for more intuitive use of the app.

4. Added calibration function.

AD Conversion and API Call Check:

- Voltage to Force Conversions
 - Function Generator Setting: Amplitude 4.00V, Vpp = 8.0V




- Result Shows that the recorded voltages are in uV, there is no ADC conversion needed. The ADC conversion could be done by the Lock-in or ActiChamp System.


- Weight > Tether Tension (force) -> Voltage -> ADC -> Digital Readings

- 300bs -> 13.078kg -> 10V -> 3036uV or 13.613kg/V

- Experimental Data with some friction, about 33lbs/V




- Lean Angle (Weight, Height, etc.) -> Tether Tension -> Voltage -> ADC -> Digital Readings -> Conversion Formula -> Actual Lean Angle



- Geometrics: R = h * sin(theta)

From Equation of Motion:

$$mg \sin(\theta) = F$$



g. Things to Check before Data Collection:

- i. Sensor values are plotted in blue and markers are plotted as orange xlines.
- ii. You want to select the appropriate channel for the accelerometer of the force sensor. The channel numbers follow the following convention:
 - 1) 32 or 64 EEG channels are channel 1 : 32 : 1 : 64.
 - 2) The specific number of AUX channels (we are using ActiChamp or LiveKamp STE, both have 8 AUX channels) are 32 : 8 or 64 + 8 => 33 : 40 : 65 : 72. The order is from AUX 1 to AUX 8.
 - 3) The LiveKamp onboard accelerometer, if enabled, will be additional 3 channels following the AUX channels, Thus 41 : 48 or 73 : 75: the order follows x, y, z.
- iv. After selecting the appropriate channel, we can check if we have the correct signal:
 - 1) If we are plotting the force sensor, try pull the tether, the readouts should change correspondingly.
 - 2) If we are plotting the accelerometer, try moving the LiveKamp or the standalone sensor.
 - 3) If we are also plotting the trigger button, click the button and see if the orange xlines, as shown in the figure above.
- v. After confirming the correct signal, we can start the output and restart the MATLAB LiveKamp Connector, and the Lab Recorder.
- h. Record: now go back to the Lab Recorder and click Start. We are ready for the experiment.
- i. Indicator Lights: the indicator light indicate whether the subject leans to a sufficient angle (specified by the Target Angle (deg)) and ready to record. Red means the angle is too large, green means the angles are too small. Green means just right. The same rule can apply for the torque by switching the Estimate type dropdown in the Lean Torque.
- j. Calibrate: Not tested yet. You should be able to calibrate the sensor using a known weight.
 - i. Connect the known weight (could be nothing) to the tethered lean-release system.
 - ii. Input the calib weight to the edit field.
 - iii. Start the stream.
 - iv. Click Calibrate. This button should be enabled.
 - v. Click Calibrate.
 - vi. An Offset will be calculated according to the difference between the known weight force and the force sensor force.
 - % Button pushed function: Calibratebutton
 - function CalibratebuttonPushed(app, event)
 % Button pushed function: Calibratebutton
 % Button pushed function: CalibratebuttonPushed(app, event)
 % The difference between the measured Force and the calibration weight force should be the offset.
 app.offset = (app.CalibWeightGdfield.Value * app.grav_acc) - app.force;
 end